## QC 623 Blood Gas•Electrolyte Control

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#### ENGLISH

**INTENDED USE**

RNA Medical® Brand QC 623 Blood Gas•Electrolyte Control is an assayed quality control material used for monitoring the performance of blood gas and electrolyte instrumentation for the analytes and analyzers listed on the Expected Values Chart.

**PRODUCT DESCRIPTION**

QC 623 is provided in three (3) levels for monitoring analyzer performance at different points within the clinical range. Level 4 is provided for elevated pO₂ monitoring. QC 623 is packaged in sealed glass ampuls, each containing 2.5 mL of solution. Ampuls are packaged thirty (30) per box.

*Active Ingredients:*

QC 623 is a buffered aqueous solution containing electrolytes (Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, Ca<sup>++</sup>). It has been equilibrated with specific levels of CO₂, O₂, and N₂. This control contains no preservatives and no human or biological materials.

**STORAGE**

The expiration date stated on the QC 623 packaging is for product stored at room temperature (up to 25 °C). Avoid freezing and temperatures greater than 30 °C.

**DIRECTIONS FOR USE**

The control should be brought to a temperature of 20-25 °C before use (see instructions regarding Expected Values). Allow at least four (4) hours for the ampuls to equilibrate to this temperature prior to testing.

For pH/blood gas values, the control should be analyzed immediately after opening. For electrolyte measurements, QC 623 is stable for up to one (1) hour after opening.

Before use, hold the ampul at the top and bottom (with forefinger and thumb) and shake for 10 seconds to mix the solution. Tap the ampul to restore the liquid to the bottom. Use gauze, tissue, gloves, or an appropriate ampul opener to protect fingers from cuts and open the ampul by snapping off the tip at the score. Immediately introduce the liquid from the ampul to the analyzer, following the instrument manufacturer's instructions for sampling a control material. Use direct aspiration, syringe transfer, or capillary mode techniques.

**EXPECTED VALUES**

The values for each control analyte on the enclosed Expected Values Chart are based on multiple determinations performed on randomly selected samples from each lot. The listing for each instrument represents the expected range and mean value of this range for ampuls that are at 25 °C when tested. (Note: pO₂ values will vary inversely by about one percent (1%) per degree Celsius that the temperature of the ampul varies from 25 °C.)

The Expected Values are provided as a guide in evaluating analyzer performance. Since instrument design and operating conditions may vary, each laboratory should establish its own expected values and control limits. The mean value established should fall within the Expected Value range shown on the chart.

**STATISTICAL SUPPORT**

RNA Medical PeerQC®, available at www.RNAMedical.com, provides monthly statistical reports for tracking and review of analyzer performance as well as lot number specific peer group data. This service is available at no charge to RNA Medical customers. Please contact RNA Medical or visit our website for information about utilizing PeerQC for this product.

**LIMITATIONS**

- QC 623 is sensitive to many instrument related factors that would affect analytical results. Because it is not a blood-based material, it may not detect certain malfunctions that would affect the testing of blood.
- This product is intended for use as a quality control material and can assist in evaluating the performance of laboratory instruments. It is not for use as a calibration standard and its use should not replace other aspects of a complete quality control program.

RNA Medical is a registered trademark and PeerQC is a registered service mark of Bionostics, Inc. The product described herein is covered by the following U.S. Patent and its foreign counterpart: 7,027,931

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LifeHealth, LLC, Roseville, MN
Medica Corporation, Bedford, MA
Siemens Healthcare Diagnostics Inc., Deerfield, IL
Via Medical, Div. of International Biomedical Ltd., Austin, TX
RNA Medical, Division of Bionostics, Inc.