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© 2010 Becton Dickinson and Company. All rights reserved. Printed in the USA.

ENGLISH

INTENDED USE

RNA Medical® Brand Hematocrit Control is a quality control material used for verifying the performance of hematocrit sensors in i-STAT cartridges.

PRODUCT DESCRIPTION

RNA Medical Hematocrit Control is provided in three (3) levels for monitoring cartridge performance at different points within the clinical range. It is packaged in sealed glass ampuls, each containing 1.7 mL of solution. Ampuls are packaged ten (10) per box.

Active Ingredients:

This product is a buffered aqueous solution containing electrolytes and non-conductive ingredients. This control contains no red cells and no human or biological materials.

STORAGE

The expiration date stated on the packaging is for product stored at temperatures of 2-25 °C. Avoid exposure to freezing and temperatures greater than 30 °C.

EXPECTED VALUES

Refer to the Abbott Point of Care Inc. website (www.abottpointofcare.com) for lot specific Value Assignment Sheets.

DIRECTIONS FOR USE

Refer to the i-STAT System Manual for use of the analyzer and cartridges.

If stored refrigerated, the control material should be equilibrated to room temperature for at least four (4) hours prior to testing.

- Before use, gently invert the ampul to mix the solution. Tap the ampul to restore the liquid to the bottom.
- Open the ampul by snapping off the tip at the score. Use gauze, tissue, gloves, or an appropriate ampul opener to protect fingers from cuts.
- Introduce the liquid from the ampul to the cartridge following the instructions for sampling a hematocrit control material outlined in the Quality Control section of the i-STAT System Manual.

***Note:** The Control pathway from the Quality Tests Menu must be used on the i-STAT 1 Analyzer.*

- Appropriately discard the ampul after use.
- Compare your result to the assigned value range. If your result is within the assigned value range, use the cartridges as needed.
- If your result is outside of the assigned value range, consult the Quality Control section of the i-STAT System Manual for troubleshooting recommendations.

***Note:** The only assigned value for this control is **hematocrit**. All other cartridge analyte results obtained with this material should be ignored.*

STATISTICAL SUPPORT

RNA Medical PeerQC™, available at www.RNAMedical.com, provides monthly statistical reports for tracking and review of cartridge performance as well as lot number specific peer group data. This service is available at no charge to RNA Medical customers. Please contact RNA Medical or visit our website for information about utilizing PeerQC for this product.

LIMITATIONS

- RNA Medical Hematocrit Control is sensitive to many instrument related factors that would affect analytical results. Because it does not contain red cells and is not a blood-based material, its ability to detect certain instrument malfunctions that would affect the testing of blood may be limited by matrix differences between the controls and blood samples.
- This product may turn slightly yellow in color over its shelf life if stored at room temperature. This coloration is normal and does not affect product performance.
- This product is intendeded for use in evaluating the performance of hematocrit sensors in i-STAT cartridges. It is not for use as a calibration standard, and its use should not replace other aspects of a complete quality control program.
- RNA Medical Hematocrit Control is intended for use on i-STAT cartridges that measure hematocrit. It is not suitable for use on analyzers that do not utilize an electrical conductivity method for hematocrit determination.

i-STAT is a registered trademark of Abbott Laboratories, Abbott Park, IL, USA. RNA Medical is a registered trademark of Bionostics, Inc., Devens, MA, USA. PeerQC is a service mark of Bionostics, Inc., Devens, MA, USA.

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Se Abbott Point of Care Inc. websitet (www.abottpointofcare.com) for lot-specifikke værditildelingsskemaer.