



January 7, 2021

Dear Valued Customer,

Our CVC123 lot 012325 now has values for the Medica EasyStat 300 Blood Gas Analyzer. Please see below for the instrument value assignments for this lot.

CVC123 Kit Lot: 012325 Exp: 2023-10

	Level 1 Lot #04123			Level 2 Lot #01533			Level 3 Lot #01632			Level 4 Lot #01733			Level 5 Lot #04224		
	Mean	Range		Mean	Range		Mean	Range		Mean	Range		Mean	Range	
Na ⁺	85.5	80.5	90.5	113.1	108.1	118.1	134.4	129.4	139.4	160.0	155.0	165.0	170.7	165.7	175.7
K ⁺	11.70	9.70	13.70	1.91	1.41	2.41	4.44	3.94	4.94	6.43	5.93	6.93	1.49	0.99	1.99
Cl ⁻	65.5	60.5	70.5	75.5	70.5	80.5	80.2	75.2	85.2	119.7	114.7	124.7	130.5	125.5	135.5
Ca ⁺⁺	2.82	2.42	3.22	1.45	1.35	1.55	1.13	1.03	1.23	0.55	0.45	0.65	0.33	0.23	0.43
pH	6.712	6.682	6.742	7.129	7.099	7.159	7.478	7.448	7.508	7.700	7.670	7.730	7.947	7.897	7.997
pO ₂	17	5	29	60	50	70	110	100	120	163	151	175	497	447	547
pCO ₂	99.6	87.6	111.6	73.6	66.6	80.6	40.6	35.6	45.6	19.3	15.3	23.3	11.0	7.0	15.0
Gluc	>500	N/A	N/A	71	62	80	192	172	212	282	255	310	<1	N/A	N/A
Lact	16.57	16.00	17.20	0.54	0.35	0.75	2.39	2.10	2.60	6.72	6.22	7.22	11.64	11.00	12.20

If you have any questions, please feel free to reach out using the contact information below.

RNA Medical, Division of Bionostics, Inc.
7 Jackson Road
Devens, MA 01434
Phone: 978-772-9070 / 800-533-6162
E-Mail: customerservice@rnamedical.com

BIOMKT-221 Rev A