

Safe-Wrap

Safe-Wrap tubes are Mylar®-wrapped glass capillary blood collection devices that have been treated with calcium-balanced lithium heparin. The Mylar contains the glass and the sample in the event of accidental breakage, meeting regulatory requirements for safety. Safe-Wrap tubes may be used for the collection of whole blood samples for blood gas, electrolyte and other testing.



Combo Blood Collection Tubes

- For use with the i-STAT® System
- Adjustable fill (65 or 95 µL cartridges)
- Plunger for easy dispensing



Blood Collection Tubes

- For use with most blood gas systems
- Markings at 55, 110 and 165 µL
- 220 µL maximum fill volume

RNA Medical also offers: blood gas/critical care controls, hemoglobin controls, calibration verification and linearity materials and hemoglobin A_{1C} controls

Please contact us for more information

RNA[®]
MEDICAL

(800) 533-6162 • +1.978.772.9070
www.RNAMedical.com

CE Labeled

