



ASSAY ASSIGNMENT APPROVAL SHEET

Part number: CVC323  
Lot number: 011601R

Level 1		Lot #	911701	Value Assignment		
		Exp Date	Dec-11			
Instrument		units	Mean	Range		
OPTI LION	Sodium	mmol/L	102	97	- 107	
	Potassium	mmol/L	2.2	1.7	- 2.7	
	Chloride	mmol/L	83	77	- 89	
	ionized Calcium	mmol/L	0.51	0.36	- 0.66	
	pH		7.72	7.66	- 7.78	

Level 2		Lot #	911801	Value Assignment		
		Exp Date	Dec-11			
Instrument		units	Mean	Range		
OPTI LION	Sodium	mmol/L	118	113	- 123	
	Potassium	mmol/L	3.0	2.5	- 3.5	
	Chloride	mmol/L	92	86	- 98	
	ionized Calcium	mmol/L	1.23	1.08	- 1.38	
	pH		7.49	7.43	- 7.55	

Level 3		Lot #	911301	Value Assignment		
		Exp Date	Dec-11			
Instrument		units	Mean	Range		
OPTI LION	Sodium	mmol/L	136	131	- 141	
	Potassium	mmol/L	4.6	4.1	- 5.1	
	Chloride	mmol/L	108	100	- 116	
	ionized Calcium	mmol/L	1.67	1.52	- 1.82	
	pH		7.29	7.23	- 7.35	

Level 4		Lot #	911401	Value Assignment		
		Exp Date	Dec-11			
Instrument		units	Mean	Range		
OPTI LION	Sodium	mmol/L	159	152	- 166	
	Potassium	mmol/L	7.1	6.6	- 7.6	
	Chloride	mmol/L	128	120	- 136	
	ionized Calcium	mmol/L	1.96	1.81	- 2.11	
	pH		7.07	7.01	- 7.13	

Level 5		Lot #	911501	Value Assignment		
		Exp Date	Dec-11			
Instrument		units	Mean	Range		
OPTI LION	Sodium	mmol/L	177	170	- 184	
	Potassium	mmol/L	9.00	8.5	- 9.5	
	Chloride	mmol/L	143	135	- 151	
	ionized Calcium	mmol/L	2.69	2.39	- 2.99	
	pH		6.89	6.83	- 6.95	

Marketing B. Harns 3/14/11  
 Logistics James O'Brien 3/14/11  
 R & D Paul T. 3/14/11  
 Quality Jan Owens 3/14/11