

RNA QC 623 Blood Gas Electrolyte Control Level 1**LOT 70217**

Analyzer	pH	PCO2 mmHg	pO2 mmHg	Ca++ mmol/l	Na+ mmol/l	K+ mmol/l	Cl- mmol/l
ECOLYTE	7.17			1.50	116	2.0	81
ECOSYSII	7.17	70	69				
COMBISYSII	7.17	70	69	1.50	116	2.0	80
MODULAR / MODULA	7.16	70	69	1.50	116	2.0	81
COMBILINE	7.17	70	69	1.50	116	2.0	80

RNA QC 623 Blood Gas Electrolyte Control Level 2**LOT 70317**

Analyzer	pH	PCO2 mmHg	pO2 mmHg	Ca++ mmol/l	Na+ mmol/l	K+ mmol/l	Cl- mmol/l
ECOLYTE	7.42			1.18	133	4.4	95
ECOSYSII	7.42	42	106				
COMBISYSII	7.42	43	106	1.18	133	4.4	95
MODULAR / MODULA	7.41	43	105	1.18	134	4.4	95
COMBILINE	7.43	43	106	1.18	133	4.4	95

RNA QC 623 Blood Gas Electrolyte Control Level 3**LOT 70417**

Analyzer	pH	PCO2 mmHg	pO2 mmHg	Ca++ mmol/l	Na+ mmol/l	K+ mmol/l	Cl- mmol/l
ECOLYTE	7.63			0.52	156	6.2	125
ECOSYSII	7.63	21	144				
COMBISYSII	7.63	21	144	0.52	156	6.2	125
MODULAR / Modular I	7.62	21	144	0.52	156	6.2	125
COMBILINE	7.64	22	144	0.52	156	6.3	125