

RNA QC 823 Range Blood Gas Electrolyte Metabolite Control Level 5**LOT 74211**

Analyzer	pH	PCO2 mmHg	pO2 mmHg	Ca++ mmol/l	Na+ mmol/l	K+ mmol/l	Cl- mmol/l	Glucose mmol/l	Lactate mmol/l
ECOLYTE	7.79			0.20	167	1.7	130		
ECOSYSII	7.79	13	430						
COMBISYSII	7.79	13	430	0.20	167	1.7	129		
MODULAR / MODULAR I	7.78	13	430	0.20	167	1.7	129	0.1	10.7
COMBILINE	7.79	13	430	0.19	167	1.7	129	0.1	10.7

BIOIVAWB-0039
Rev. A